

LID 574 Toolbox set up with BLE Base Unit





Welcome to the GR252, LID574, LID674 settings

This wizard will guide you through the
GR252, LID574, LID674 settings.

Using this wizard you can edit the settings
of your handheld scanner.

To continue, click Next.



< Back

Next >

Cancel

Select transponders to read

Please select the transponder types the reader should read



Current firmware version: LID574V1.1.1CR-B

Trovan unique

Trovan flex 3 blocks

PSK1

Trovan flex 5 blocks

PSK2

Trovan flex 7 blocks

To increase reading performance only select required transponder types.



< Back

Next >

Cancel

General settings

Select general settings



Use custom names

Show only custom name

Read only new ID

Show last code

Shutdown time after read (s):

Trigger time (s):

Buzzer time after read (s):

Haptic feedback time after read (s):

Buzzer time custom name found (s):

Haptic feedback time custom name found (s):



< Back

Next >

Cancel

→ GR252, LID574, LID674 settings writer V105



Save settings

Save tag code to reader memory options



Save codes

Memory size (KB) 2048

with date/time

Up to number of Trovan codes without date/time: 95325

Up to number of Trovan codes with date/time: 61680

Memory used: 0 %



< Back

Next >

Cancel

GR252, LID574, LID674 settings writer V105



Communication settings

Set communication options



Send id by USB

Send ID by USB

USB Comport

USB Keyboard

USB CDC format

<ID><CR>



< Back

Next >

Cancel

Communication settings: BLE

Wireless connection settings



Send id by Bluetooth

- Send ID by Bluetooth
- Keep connected
- Always connected
- Bootup connect, after sent disconnect
- Connect when tag found, after sent disconnect

Ble format

<ID>

ID

Tag type HEX Tag type text Custom name

<CR> <LF> / <ENTER> <TAB>



< Back Next > Cancel

→ GR252, LID574, LID674 settings writer V105



Completed
Ready for use



The reader is set up. Please disconnect the device.
Please disconnect the reader.



< Back

Finish


Cancel

Connecting the BLE-Base to your Computer


1. Plug the BLE-Base USB cable to your computer
2. Press the RED Button on the BLE-Base
3. Double click on the LID574- The display will state connecting, then connected
4. Please see additional documentation on connecting the BLE-Base in (BLE BASE Manual)

Connect to Device Manager on your computer


Double click on USB Serial Port (COM ?)



1. If the port number does not show COM3 (for Provantis) the click on Advance



2. Change the Port number to COM 3



Advanced Settings for COM3

? X

COM Port Number: COM3

OK

USB Transfer Sizes

Select lower settings to correct performance problems at low baud rates.

Select higher settings for faster performance.

Receive (Bytes): 4096

Transmit (Bytes): 4096

Cancel

Defaults

BM Options

Select lower settings to correct response problems.

Latency Timer (msec): 16

Timeouts

Minimum Read Timeout (msec): 0

Minimum Write Timeout (msec): 0

Miscellaneous Options

Serial Enumerator

Serial Printer

Cancel If Power Off

Event On Surprise Removal


Set RTS On Close

Disable Modem Ctrl At Startup


Enable Selective Suspend

Selective Suspend Idle Timeout (secs): 5

3. Note the changes stating COM 3 In Device Manager



Testing the Read Using the LID574 and Bluetooth transmission to the BLE-Base



Connected to: COM3 | USB Serial Port , Comtest RS232

Comport: COM3 | USB Serial Port

Connect Disconnect

Baudrate: 110 300 600 1200 2400 4800 9600 14400 19200 38400 56000 57600 115200

Parity: none odd even mark space

Handshaking: none RTS/CTS XON/XOFF RTS/CTS + XON/XOFF

Databits: 5 6 7 8

Receive CLEAR HEX Add date/time Receive end of line: <CR>

```
00071F8DE2<CR>
0006CE97E9<CR>
```

Transmit CLEAR SET DTR CLR DTR SET RTS CLR RTS Transmit end of line: <CR>

Connected to: COM3 | USB Serial Port Transmitted:

Images

